



Recombinant Asparagine--tRNA ligase (asnS)

Product Code	CSB-BP349169EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P58694
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	SVVPVADVL QGRVAVDSEV TVRGWVRTRR DSKAGISFLA VYDGSCFDPV QAVINNSLPN YNEDVLRLLT GCSVIVTGKV VASPGQGQQF EIQTSKVEVA GWVEDPDTYP MAAKRHSIEY LREVAHLRPR TNLIGAVARV RHTLAQALHR FFNEQGFFWV STPLITASDT EGAGEMFRVS TLDLENLPRN DQGKVDFFDKD FFGKESFLTV SGQLNGETYA CALSKIYTFG PTFRAENSNT SRHLAEFWML EPEVAFANLN DIAGLAEAML KYVFKAVLEE RADDMKFFAE RVDKDAVSRL ERFIEADFAQ VDYTEAVTIL ENCGRKFNVP VYWGVDLSSE HERYLAEHF KAPVVVKNYP KDIKAFYMRL NEDGKTVAAM DVLAPGIGEI IGGSQREERL DVLDERMLEM GLNKEDYWWY RDLRRYGTVP HSGFGLGFER LIAYVTGVQN VRDVIPFPRT PRNASF
Source	Baculovirus
Target Names	asnS
Protein Names	Recommended name: Asparagine--tRNA ligase EC= 6.1.1.22 Alternative name(s): Asparaginyl-tRNA synthetase Short name= AsnRS
Expression Region	2-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.