



# Recombinant Human Tubby-related protein 2 (TULP2)

<b>Product Code</b>	CSB-MP025347HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O00295
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSQDNNTLMR DILGHELAAM RLQKLEQQRR LFEKKQRQKR QELLMVQANP DASPWLWRSC LREERLLGDR GLGNPFLRKK VSEAHLPSCI HSALGTVSCG GDGRGERGLP TPRTEAVFRN LGLQSPFLSW LPDNSDAELE EVSVENGVS PPPFKQSPRI RRGWQAHQR PGTRAEGESD SQDMGDAHKS PNMGPNPGMD GDCVYENLAF QKEEDLEKKR EASESTGTNS SAAHNEELSK ALKGEGGTDS DHMRHEASLA IRSPCPGLEE DMEAYVLRPA LPGTMMQCYL TRDKHGVKDG LFPLYLYLE TSDSLQRFL AGRKRRRSKT SNYLISLPT HLSRDGDNFV GKVRSNVFST KFTIFDNGVN PDREHLTRNT ARIRQELGAV CYEPNVLYL GPRKMTVILP GTNSQNQRIN VQPLNEQESL LSRVQRGDKQ GLLLLHNKTP SWDKENGVYT LNFHGRVTRA SVKNFQIVDP KHQEHLVLQF GRVGPDTFTM DFCFPFSPLQ AFSICLSSFN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TULP2
<b>Protein Names</b>	Recommended name: Tubby-related protein 2 Alternative name(s): Cancer/testis antigen 65 Short name= CT65 Tubby-like protein 2
<b>Expression Region</b>	1-520
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	TULP2 is a member of a family of tubby-like genes (TULPs) that encode proteins of unknown function. Members of this family have been identified in plants, vertebrates, and invertebrates. The TULP proteins share a conserved C-terminal region of approximately 200 amino acid residues.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.