



# Recombinant Cricetulus longicaudatus Guanine nucleotide-binding protein G (o) subunit alpha

<b>Product Code</b>	CSB-EP349372DXV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P59216
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Cricetulus longicaudatus (Long-tailed dwarf hamster) (Chinese hamster)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GCTLSAEER AALERSKAIE KNLKEDGISA AKDVKLLLLG AGESGKSTIV KQMKIIHEDG FSGEDVKQYK PVVYSNTIQS LAAIVRAMDT LGVEYGDKER KADSKMVCDV VSRMEDTEPF SAELLSAMMR LWGDSGIQEC FNRSREYQLN DSAKYYLDSL DRIGAADYQP TEQDILRTRV KTTGIVETHF TFKNLHFRLF DVGGQRSEK KWIHCFEDVT AIFCVALSG YDQVLHEDET TNRMHESLML FDSICNNKFF IDTSIILFLN KKDLFGEKIK KSPLTICFPE YPGSNTYEDA AAYIQTFES KNRSPNKEY CHMTCATDTN NIQVVFDAVT DIIANNLRG CGLY
<b>Source</b>	E.coli
<b>Target Names</b>	GNAO1
<b>Protein Names</b>	Recommended name: Guanine nucleotide-binding protein G(o) subunit alpha
<b>Expression Region</b>	2-354
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.