



Recombinant Pig Fatty acid-binding protein, heart (FABP3)

Product Code	CSB-BP007943PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02772
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	VDAFAGTWK LVDSKNFDDY MKSIGVGFAT RQVANMTKPT TIEVNGDTI IIKTQSTFKS TEISFKLGVE FDETTADDRK VKSIVTLDGG KLVHLQKWNG QETTLVRELV DGKLILTLTH GSAVCTR TYE KEA
Source	Baculovirus
Target Names	FABP3
Protein Names	Recommended name: Fatty acid-binding protein, heart Alternative name(s): Fatty acid-binding protein 3 Heart-type fatty acid-binding protein Short name= H- FABP
Expression Region	2-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The intracellular fatty acid-binding proteins (FABPs) belongs to a multigene family. FABPs are divided into at least three distinct types, namely the hepatic- intestinal- and cardiac-type. They form 14-15 kDa proteins and are thought to participate in the uptake, intracellular metabolism and/or transport of long-chain fatty acids. They may also be responsible in the modulation of cell growth and proliferation. Fatty acid-binding protein 3 gene contains four exons and its function is to arrest growth of mammary epithelial cells. This gene is a candidate tumor suppressor gene for human breast cancer.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.