



Recombinant Buchnera aphidicola subsp. Acyrthosiphon pisum SsrA-binding protein (smpB)

Product Code	CSB-BP349648BWZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P57342
Product Type	Recombinant Protein
Immunogen Species	Buchnera aphidicola subsp. Acyrthosiphon pisum (strain APS) (Acyrthosiphon pisum symbiotic bacterium)
Purity	>85% (SDS-PAGE)
Sequence	MLQKKKYQKK SSKIIIINKKA YYNYFIEKVF QSGLVLEGWE IKSIRSGKVN ISESYIINDR NEMYLCNCLI EPLQMSSNRF SCDPTRKKKL LLHKNEIDFL SLKKKNTGYT MVSLSLFWKK SWCKLEFGLA KGKTAQDKRI NLKKREWEQE KLKILKKTKE TY
Source	Baculovirus
Target Names	smpB
Protein Names	Recommended name: SsrA-binding protein
Expression Region	1-162
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.