



Recombinant Mouse Hepatocyte nuclear factor 6 (Onecut1)

Product Code	CSB-BP016336MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O08755
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNAQLTMEAI GELHGV SHEP VPAPADLLGG SPHARSSVGH RGS HLPPAHP RSMGMASLLD GGSGGSDYHH HHRAP EHS LA GPLHPTMTMA CETPPGMSMP TTYTTLPLQ PLPISTVSD KFP HHHHHHH HHHHP HHHQR LAGNVSGSFT LMRDERGLAS MNNLYTPYHK DVAGMGQSLS PLSG SGLGSI HNSQQGLPHY AHPGAAMPTD KMLTPNGFEA HHPAMLGRHG EQHLTPTSAG MVPINGLPPH HPHAHLNAQG HGQLLG TARE PNPSVTGAQV SNGSNSGQME EINTKEVAQR ITTELKRYSI PQAIFAQRVL CRSQGTLS DL LRNPKPWSKL KSGRETFRRM WKWLQEPEFQ RMSALRLAAC KRKEQEHGKD RGNTPKK PRL VFTDVQRRTL HAIFKENKRP SKELQITISQ QLGLELSTVS NFFMNARRRS LDKWQDEGGS NSGSSSSSSS TCTKA
Source	Baculovirus
Target Names	Onecut1
Protein Names	Recommended name: Hepatocyte nuclear factor 6 Short name= HNF-6 Alternative name(s): One cut domain family member 1 One cut homeobox 1
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.