



# Recombinant Mouse Netrin-1 (Ntn1)

<b>Product Code</b>	CSB-BP016127MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O09118
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           GPGLSM FAGQAAQPDP CSDENGHPRR CIPDFVNAAF GKDVRVSSTC            GRPPARYCVV SERGEERLRS CHLCNSSDPK KAHPPAFLTD LNNPHNLTCW            QSENYLQFPH NVTLTSLGK KFEVTVVSLQ FCSPRPESMA IYKSMGYGRT            WVPFQFYSTQ CRKMYNRPHR APITKQNEQE AVCTDSHTDM RPLSGGLIAF            STLDGRPSAH DFDNSPVLQD WVTATDIRVA FSRLHTFGDE NEDDSELARD            SYYYAVSDLQ VGGRCCKNGH AARCVRDRDD SLVCD CRHNT            AGPECDRCKP FHYDRPWQRA TAREANECVA CNCNLHARRC            RFNMELYKLS GRKSGGVCLN CRHNTAGRHC HYCKEGFYRD MGKPITHRKA            CKACDCHPVG AAGKTCNQTT GQCPCKDGVG GITCNRCAKG YQQSRSPIAP            CIKIPVAPPT TAASSVEEPE DCDSYCKASK GKLKMNMKKY CRKDYAVQIH            ILKADKAGDW WKFTVNIISV YKQGTSRIRR GDQSLWIRSR DIACKCPKIK            PLKKYLLLG N AEDSPDQSGI VADKSSSLVIQ WRDTWARRLR KFQQREKKGK            CKKA         </p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	Ntn1
<b>Protein Names</b>	Recommended name: Netrin-1
<b>Expression Region</b>	25-604
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Netrin is included in a family of laminin-related secreted proteins. The function of this gene has not yet been defined; however, netrin is thought to be involved in axon guidance and cell migration during development. Mutations and loss of expression of netrin suggest that variation in netrin may be involved in cancer development.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.