



Recombinant Helicobacter pylori Transcription termination factor Rho (rho)

Product Code	CSB-EP349984HUV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56466
Product Type	Recombinant Protein
Immunogen Species	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Purity	>85% (SDS-PAGE)
Sequence	MNENAPTHKS SHKVKTHTPV SGYHIEDLRT YPTEKLLEIA NKLKVENPQE FKRQDLMFEI LKTQVTQGGY ILFTGILEIM PDGYGFLRGF DGSFSDGHND TYVSPSQIRR FALRNGDIVT GQVRSPKDQE KYYALLKIEA INYSPSDEIR NRPLFDNLTP LFPDEQIKLE YEPTKVTGRM LDLFSPVGKG QRALIVAPPR TGKTELMKEL AQQITSNHPE VELIILLVDE RPEEVTDMMQR SVKGQVFSST FDLPANNHIR IAELVLERAK RRVEMGKDVV VLLDSITRLA RAYNAVTPSS GKVLSSGGVDA NALHRPKRFF GAARNIEEGG SLTIIATALI ETGSRMDEVI FEEFKGTGNS EIVLARNIAD RRIYPAFDIL KSGTRKDNIL LGKDRLTKVW VLRNVMQQMD DIEALSFVYS KMQQTKDNEE FLNLMNEK
Source	E.coli
Target Names	rho
Protein Names	Recommended name: Transcription termination factor Rho EC= 3.6.4.- Alternative name(s): ATP-dependent helicase Rho
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.