



Recombinant Human TRAF-type zinc finger domain-containing protein 1 (TRAFD1)

Product Code	CSB-BP024156HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14545
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AEFLDDQET RLCDNCKKEI PVFNFTIHEI HCQRNIGMCP TCKEPFPKSD METHMAAEHC QVTCKCNKKL EKRLKKHEE TECPLRLAVC QHCDLELSIL KLKEHEDYCG ARTELCGNCG RNVLVKDLKT HPEVCGREGE EKRNEVAIPP NAYDESWGQD GIWIASQLLR QIEALDPPMR LPRRPLRAFE SDVFHNRTTN QRNITAQVSI QNNLFEEQER QERNRGQQPP KEGGEESANL DFMLALSQN EGQASSVAEQ DFWRAVCEAD QSHGGPRSLS DIKGADEIM LPCEFCEELY PEELLIDHQT SCNPSRALPS LNTGSSSPRG VEEDVIFQN FLQQAASNQL DSLMGLSNSH PVEESIIIPC EFCGVQLEEE VLFHHQDQCD QRPATATNHV TEGIPRLDSQ PQETSPELPR RRVRHQGDLS SGYLDLTKQE TANGPTSCLP PSRPINMTA TYNQLSRSTS GPRPGCQPSS PCVPKLSNSD SQDIQGRNRD SQNGAIAPGH VSVIRPPQNL YPENIVPSFS PGPSGRYGAS GRSEGGRNSR VTPAAANYRS RTAKAKPSKQ QGAGDAEEEE EE
Source	Baculovirus
Target Names	TRAFD1
Protein Names	Recommended name: TRAF-type zinc finger domain-containing protein 1 Alternative name(s): Protein FLN29
Expression Region	2-582
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.