



# Recombinant Human TRAF-type zinc finger domain-containing protein 1 (TRAFD1)

<b>Product Code</b>	CSB-MP024156HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O14545
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AEFLDDQET RLCDNCKKEI PVFNFTIHEI HCQRNIGMCP TCKEFPKSD METHMAAHC QVTCKCNKKL EKRLKKHEE TECPLRLAVC QHCDLELSIL KLKEHEDYCG ARTELCGNCG RNVLVKDLKT HPEVCGREGE EKRNEVAIPP NAYDESWGQD GIWIASQLLR QIEALDPPMR LPRRPLRAFE SDVFHNRTTN QRNITAQVSI QNNLFEEQER QERNRGQQPP KEGGEESANL DFMLALSQN EGQASSVAEQ DFWRAVCEAD QSHGGPRSLS DIKGADEIM LPCEFCEELY PEELLIDHQT SCNPSRALPS LNTGSSSPRG VEEPDIQFN FLQQAASNQL DSLMGLSNH PVEESIIIPC EFCGVQLEEE VLFHHQDQCD QRPATATNHV TEGIPRLDSQ PQETSPQLPR RRVRHQGDLS SGYLDLTKQE TANGPTSCLP PSRPINMTA TYNQLSRSTS GPRPGCQPS PCVPKLSNSD SQDIQGRNRD SQNGAIAPGH VSVIRPPQNL YPENIVPSFS PGPSGRYGAS GRSEGGRNSR VTPAAANYRS RTAKAKPSKQ QGAGDAEEEE EE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TRAFD1
<b>Protein Names</b>	Recommended name: TRAF-type zinc finger domain-containing protein 1 Alternative name(s): Protein FLN29
<b>Expression Region</b>	2-582
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.