



Recombinant Probable transcriptional regulatory protein sgaR (sgaR)

Product Code	CSB-MP349994HYC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56644
Product Type	Recombinant Protein
Immunogen Species	Hyphomicrobium methylovorum
Purity	>85% (SDS-PAGE)
Sequence	MVAPMASKIE VSIVDQNPVV RTGLETLLVR DGRFSVSGIY ASGEDLLAAL QTKPVEIVIV GWSLPDMTGG DVLTRVRKEK WRTRIIIYTG ERSSEVLRTA IKSGAWGFVA KTEEPQVLE AVVSVARGRL SLPYVDIDLL NHDPLESLTA RERELLAALA NGWTNLQIAA RTGISRNTVK YHLKNLYDKL GVSNRAMAVA LHVSINRNEH R
Source	Mammalian cell
Target Names	sgaR
Protein Names	Recommended name: Probable transcriptional regulatory protein sgaR
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.