



# Recombinant Human DNA-binding protein RFXANK (RFXANK)

<b>Product Code</b>	CSB-EP019617HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O14593
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MELTQPAEDL IQTQQTPASE LGDPEDPGEE AADGSDTVVL SLFPCTPEPV NPEPDASVSS PQAGSSLKHS TTLTNRQRGN EVSALPATLD SLSIHQLAAQ GELDQLKEHL RKGDNLVNKP DERGFTPLIW ASAFGEIETV RFLLEWGADP HILAKERESA LSLASTGGYT DIVGLLLERD VDINIYDWNG GTPLLYAVRG NHVKCVEALL ARGADLTTEA DSGYTPMDLA VALGYRKVQQ VIENHILKLF QSNLVPADPE
<b>Source</b>	E.coli
<b>Target Names</b>	RFXANK
<b>Protein Names</b>	Recommended name: DNA-binding protein RFXANK Alternative name(s): Ankyrin repeat family A protein 1 Regulatory factor X subunit B Short name= RFX-B Regulatory factor X-associated ankyrin-containing protein
<b>Expression Region</b>	1-260
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Major histocompatibility (MHC) class II molecules are transmembrane proteins that have a central role in development and control of the immune system. This protein, along with regulatory factor X-associated protein and regulatory factor-5, forms a complex that binds to the X box motif of certain MHC class II gene promoters and activates their transcription. Once bound to the promoter, this complex associates with the non-DNA-binding factor MHC class II transactivator, which controls the cell type specificity and inducibility of MHC class II gene expression. This protein contains ankyrin repeats involved in protein-protein interactions. Mutations in this gene have been linked to bare lymphocyte syndrome type II, complementation group B. Two transcript variants encoding different isoforms have been described for this gene, with only one isoform showing activation activity.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.