



# Recombinant Mouse La-related protein 7 (Larp7)

<b>Product Code</b>	CSB-BP012756MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q05CL8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	METENQKTME ESTKRKEEKK KRSRVKQVLA DIAKQVDFWF GDANLHKDKF LREQIEKSRD GYVDISLLVS FNKMKLLTTD GKLIARALKS SSVVELDLEG TRIRRKPLG ERPKDEEERT VYVELLPKNV THSWIERVFG KCGNVVYISI PHYKSTGDPK GFAFVEFETK EQAAKAIEFL NNPPEEAPRK PGIFPKTVKN KPIPSLRVAE EKKKKKKKKG RIKKEESVQA KESAVDSSSS GVCKATKRPR TASEGSEAET PEAPKQPAKK KKKRDKVEAS SLPEARAGKR ERCSAEDEDC LPPRPKAKKR AQKDGVGQAA SEVSKESRDL EFCSTEEKE TDRKGDLSLK VKRKHKKKHK ERHKMGEEVI PLRVLSKTEW MDLKKEYLAL QKASMALKK TISQIKLESE METDCKAPTA GSGQECSTQE KVS AQGPQFV TGVIVKIVSG EPLPGRKQVK DILATISEVV YIDLLEGDTE CHARFKTPED AQAVMNAQTE IRKKHSWNLE VLSGDHEQRY WQKILVDRQA KLNQPREKKR GTEKLITKAE KIRLAKTQQA SQHIRFSEYD
<b>Source</b>	Baculovirus
<b>Target Names</b>	Larp7
<b>Protein Names</b>	Recommended name: La-related protein 7 Alternative name(s): La ribonucleoprotein domain family member 7
<b>Expression Region</b>	1-570
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.