



Recombinant Mouse La-related protein 7 (Larp7)

Product Code	CSB-MP012756MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05CL8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	METENQKTME ESTKRKEEKK KRSRVKQVLA DIAKQVDFWF GDANLHKDKF LREQIEKSRD GYVDISLLVS FNKMKKLTDD GKLIAALKS SSVVELDLEG TRIRRKPLG ERPKDEEERT VYVELLPKNV THSWIERVFG KCGNVVYISI PHYKSTGDPK GFAFVEFETK EQAAKAIEFL NNPPEEAPRK PGIFPKTVKN KPIPSLRVAE EKKKKKKKKG RIKKEESVQA KESAVDSSSS GVCKATKRPR TASEGSEAET PEAPKQPAKK KKKRDKVEAS SLPEARAGKR ERCSAEDEDC LPPRPKAKKR AQKDGVGQAA SEVSKESRDL EFCSTEEEEK TDRKGDLSLK VKRKHKKKHK ERHKMGEEVI PLRVLSKTEW MDLKKEYLAL QKASMALKK TISQIKLESE METDCKAPTA GSGQECSTQE KVSAAQGPQFV TGVIVKIVSG EPLPGRKQVK DILATISEVV YIDLLEGDTE CHARFKTPED AQAVMNAQTE IRKKHSWNLE VLSGDHEQRY WQKILVDRQA KLNQPREKKR GTEKLITKAE KIRLAKTQQA SQHIRFSEYD
Source	Mammalian cell
Target Names	Larp7
Protein Names	Recommended name: La-related protein 7 Alternative name(s): La ribonucleoprotein domain family member 7
Expression Region	1-570
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.