



# Recombinant Human Zinc finger protein 195 (ZNF195)

<b>Product Code</b>	CSB-BP026570HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O14628
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTLLTFRDVA IEFSLLEEWKC LDLAQQNLYR DVMLENYRNL FSVGLTVCKP GLITCLEQRK EPWNVKRQEA ADGHPENMGFH HATQACLELL GSSDLPASAS QSAGITGVNH RAQPGLNVSV DKFTALCSPG VLQTVKWFLE FRCIFSLAMS SHFTQDLLPE QGIQDAFPKR ILRGYGNCGL DNLYLRKDWE SLDECKLQKD YNGLNQCSST THSKIFQYNK YVKIFDNFSN LHRNINISNTG EKPFKCECG KSFQMLSFLT EHQKIHTGKK FQKCGECGKT FIQCSTFTEP ENIDTGEKPY KCQECNNVIK TCSVLTKNRI YAGGEHYRCE EFGKVFNQCS HLTEHEHGTE EKPCYEECS SVFISCSLS NQQMILAGEK LSKCETWYKG FNHSPNPSKH QRNEIGGKPF KCEECDSIFK WFSDLTKHHR IHTGEKPYKC DECGKAYTQS SHLSEHRRIH TGEKPYQCEE CGKVFRTCSS LSNHKRTHSE EKPYTCEECG NIFKQLSDLT KHKKTHHTGEK PYKCDECGKN FTQSSNLIVH KRIHTGEKPY KCEECGRVFM WFSDITKHKK THTGEKPYKC DECGKNFTQS SNLIVHKRIH TGEKPYKCEK CGKAFTQFSH LTVHESIHT
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF195
<b>Protein Names</b>	Recommended name: Zinc finger protein 195
<b>Expression Region</b>	1-629
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.