



Recombinant Human Zinc finger protein 195 (ZNF195)

Product Code	CSB-MP026570HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14628
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTLLTFRDVA IEFSLLEEWKC LDLAQQNLYR DVMLENYRNL FSVGLTVCKP GLITCLEQRK EPWNVKRQEA ADGHPENMGFH HATQACLELL GSSDLPASAS QSAGITGVNH RAQPGLNVSV DKFTALCSPG VLQTVKWFLE FRCIFSLAMS SHFTQDLLPE QGIQDAFPKR ILRGYGNCGL DNLYLRKDWE SLDECKLQKD YNGLNQCST TSHKIFQYNK YVKIFDNFSN LHRNINISNTG EKPFKCECG KSFQMLSFLT EHQKIHTGKK FQKCGECGKT FIQCSTFTEP ENIDTGEKPY KCQECNNVIK TCSVLTKNRI YAGGEHYRCE EFGKVFNQCS HLTEHEHGTE EKPCYEECS SVFISCSLS NQQMILAGEK LSKCETWYKG FNHSPNPSKH QRNEIGGKPF KCEECDSIFK WFSDLTKHHR IHTGEKPYKC DECGKAYTQS SHLSEHRRIH TGEKPYQCEE CGKVFRTCSS LSNHKRTHSE EKPYTCEECG NIFKQLSDLT KHKKTHHTGEK PYKCDECGKN FTQSSNLIVH KRIHTGEKPY KCEECGRVFM WFSDITKHKK THTGEKPYKC DECGKNFTQS SNLIVHKRIH TGEKPYKCEK CGKAFTQFSH LTVHESIHT</p>
Source	Mammalian cell
Target Names	ZNF195
Protein Names	Recommended name: Zinc finger protein 195
Expression Region	1-629
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.