



Recombinant Human Ectoderm-neural cortex protein 1 (ENC1)

Product Code	CSB-MP007663HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14682
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSVSVHENRK SRASSGSINI YLFHKSSYAD SVLTHLNLLR QQRLFTDVLL HAGNRTFPCH RAVLAACSR Y FEAMFSGGLK ESQDSEVNFD NSIHPEVLEL LLDYAYSSRV IINEENAESL LEAGDMLEFQ DIRDACAFL EKNLHPTNCL GMLLSDAHQ CTKLYELSWR MCLSNFQ TIR KNEDFLQLPQ DMVVQLLSSE ELETEDERLV YESAINWISY DLKKRYCYLP ELLQTVRLAL LPAIYLMENV AMEELITKQR KSKEIVEEAI RCKLKILQND GVVTS LCARP RKTGHALFLL GGQTFMCDKL YLVDQKAKEI IPKADIPSPR KEFSACAIGC KVIYTGGRGS ENGVSKDVWV YDTLHEEWSK AAPMLVARFG HGSAELKHCL YVVGHTAAT GCLPASPSVS LKQVEHYDPT INKWTMVAPL REGVSNAAVV SAKLKLFAFG GTSVSHDKLP KVQCYDQCEN RWTVPATCPQ PWR YTAAAVL GNQIFIMGGD TEFSACSAYK FNSETYQWTK VGDVTAKRMS CHAVASGNKL YVVG GYFGIQ RCKTLDCYDP TLDVWNSITT VPYSLIPTAF VSTWKHLPS
Source	Mammalian cell
Target Names	ENC1
Protein Names	Recommended name: Ectoderm-neural cortex protein 1 Short name= ENC-1 Alternative name(s): Kelch-like protein 37 Nuclear matrix protein NRP/B p53- induced gene 10 protein
Expression Region	1-589
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.