



Recombinant Human Integrin beta-1-binding protein 1 (ITGB1BP1)

Product Code	CSB-YP011881HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14713
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFRKGGKKRHS SSSSQSSEIS TKSKSVDSSL GGLSRSTVA SLDTDSTKSS GQSNNSDTC AEFRIKYVGA IEKLLSEGGK GLEGPLDLIN YIDVAQQDGK LPFVPEEEF IMGVSKYGIK VSTSDQYDVL HRHALYLIIR MVCYDDGLGA GKSLALKTT DASNEEYSLW VYQCNSLEQA QAICKVLSTA FDSVLTSEKP
Source	Yeast
Target Names	ITGB1BP1
Protein Names	Recommended name: Integrin beta-1-binding protein 1 Alternative name(s): Integrin cytoplasmic domain-associated protein 1 Short name= ICAP-1
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The cytoplasmic domains of integrins are essential for cell adhesion. This protein binds to the beta1 integrin cytoplasmic domain. The interaction between this protein and beta1 integrin is highly specific. Two isoforms of this protein are derived from alternatively spliced transcripts. The shorter form of this protein does not interact with the beta1 integrin cytoplasmic domain. The longer form is a phosphoprotein and the extent of its phosphorylation is regulated by the cell-matrix interaction, suggesting an important role of this protein during integrin-dependent cell adhesion.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.