



Recombinant Human Coiled-coil domain-containing protein 67 (CCDC67)

Product Code	CSB-MP004729HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05D60
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MENQAHTMG TSPCEAELQE LMEQIDIMVS NKKMDWERKM RALETRDLR DQELANAQTC LDQKGQEVGL LRQKLDSEK CNLAMTQNYE GQLQSLKAQF SKLTNNFEKL RLHQMKQNKV PRKELPHLKE EIPFELSNLN QKLEEFRAKS REWDKQEILY QTHLISLDAQ QKLLSEKCNQ FQKQAQSYQT QLNGKKQCLE DSSSEIPRLI CDPDPNCEIN ERDEFIIEKL KSAVNEIALS RNKLQDENQK LLQELKMYQR QCQAMEAGLS EVKSELQSRD DLLRIIEMER LQLHRELLKI GECQNAQGNK TRLESSYLPS IKEPERKIKE LFSVMQDQPN HEKELNKIRS QLQQVEEYHN SEQERMNEI SDLTEELHQQ EITIAVTCK AALLEKQLKM ELEIKEKMLA KQKVSDMKYK AVRTENTHLK GMMGDLDPGE YMSMDFTNRE QSRHTSINKL QYENERLRND LAKLHVNGKS TWTNQNTYEE TGRYAYQSQI KVEQNEERLS HDCEPNRSTM PPLPPSTFQA KEMTSPLVSD DDVFPLSPPD MSFPASLAAQ HFLLEEEKRA KELEKLLNTH IDELQRHTEF TLNKYSKLLKQ NRHI
Source	Mammalian cell
Target Names	DEUP1
Protein Names	Recommended name: Coiled-coil domain-containing protein 67
Expression Region	1-604
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.