



Recombinant Escherichia coli Ditrans,polycis-undecaprenyl-diphosphate synthase ((2E,6E)-farnesyl-diphosphate specific) (ispU)

Product Code	CSB-YP350241ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P60472
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MMLSATQPLS EKLPAHGCRH VAIIMDGNGR WAKKQGKIRA FGHKAGAKSV RRAVSFAANN GIEALTLYAF SSENWNRPAQ EVSALMELFV WALDSEVKSL HRHNVRLRII GDTSRFNSRL QERIRKSEAL TAGNTGLTLN IAANYGGRWD IVQGVRLAE KVQQGNLQPD QIDEEMLNQH VCMHELAPVD LVIRTGGEHR ISNFWLWQIA YAELYFTDVL WPDFDEQDFE GALNAFANRE RRFGGTEPGD ETA
Source	Yeast
Target Names	ispU
Protein Names	Recommended name: Undecaprenyl pyrophosphate synthase Short name= UPP synthase EC= 2.5.1.31Alternative name(s): Di-trans,poly-cis-decaprenylcistransferase Ditrans,polycis-undecaprenyl-diphosphate synthase Undecaprenyl diphosphate s
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.