



# Recombinant Human FYN-binding protein (FYB)

<b>Product Code</b>	CSB-BP009099HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O15117
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MAKYNTGGNP TEDVSVNSRP FRVTGPNSSS GIQARKNLFN NQGNASPPAG            PSNVPKFGSP KPPVAVKPSS EEKPDKEPKP PFLKPTGAGQ RFGTPASLTT            RDPEAKVGFL KPVGPKPINL PKEDSKPTFP WPPGNKPSLH SVNQDHDLPK            LGPKSGPTTP TSENEQKQAF PKLTGVKGF MSASQDLEPK PLFPKPAFGQ            KPPLSTENSH EDESPMKNVS SSKGSPAPLG VRSKSGPLKP ARENSENKDH            AGEISLPPF GVVLPKAASR GPGLSKNGE EKEDRKIDA AKNTFQSKIN            QEELASGTPP ARFPKAPSKL TVGGPWGQSQ EKEKGDKNSA TPKQKPLPPL            FTLGPPPPKP NRPPNVDLTK FHKTSNGNST SKGQTSYSTT SLPPPPPSHP            ASQPPLPASH PSQPPVPSLP PRNIKPPFDL KSPVNEDNQD GVTHSDGAGN            LDEEQDSEGE TYEDIEASKE REKKREKEEK KRLELEKKEQ KEKEKKEQEI            KKKFKLTGPI QVIHLAKACC DVKGGKNELS FKQGEQIEI RITDNPEGKW            LGRTARGSYG YIKTTAVEID YDSLKLLKDS LGAPSRPIED DQEVYDDVAE            QDDISSHSQS GSGGIFPPPP DDDIYDGIEE EDADDGFPAP PKQLDMGDEV            YDDVDTSDFP VSSAEMSQGT NVGKAKTEEK DLKLLKKQEK EEKDFRKKFK            YDGEIRVLYS TKVTTITSK KWGTRDLQVK PGESLEVIQT TDDTKVLCRN            EEGKYGYVLR SYLADNDGEI YDDIADGCIY DND</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	FYB1
<b>Protein Names</b>	<p>Recommended name: FYN-binding protein            Alternative name(s): Adhesion and degranulation promoting adaptor protein            Short name= ADAP FYB-120/130            Short name= p120/p130 FYN-T-binding protein SLAP-130 SLP-76-associated phosph</p>
<b>Expression Region</b>	1-783
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.