



Recombinant Arabidopsis thaliana Dehydration-responsive element-binding protein 2G (DREB2G)

Product Code	CSB-EP350389DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61827
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEEEQPPAKK RNMGRSRKGC MKGKGGPEN A TCTFRGVRQR TWGKWVAEIR EPNRCTRLWL GTFNTSVEAA MAYDEAAKKL YGHEAKLNLV HPQQQQQVVV NRNLSFSGHG SGWAYNKKL DMVHGLDLGL GQASCSRGSC SERSSFLLQED DDHSHNRCSS SSGSNLCWLL PKQSDSQDQE TVNATTSYGG EGGGGSTLTF STNLKPKNLM SQNYGLYNGA WSRFLVGQEK KTEHDVSSSC GSSDNKESML VPSCGGERMH RPELEERTGY LEMDDLLEID DLGLLIGKNG DFKNWCCEEF QHPWNWF
Source	E.coli
Target Names	DREB2G
Protein Names	Recommended name: Dehydration-responsive element-binding protein 2G Short name= Protein DREB2G
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.