



Recombinant Human CTD small phosphatase-like protein (CTDSPL)

Product Code	CSB-BP006144HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15194
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDGPAITQV TNPKEDEGRL PGAGEKASQC NVSLKKQRSR SILSSFFCCF RDYNVEAPPP SSPSVLPPLV EENGGLQKGD QRQVIPIPS PAKYLLPEVT VLDYGKKCVV IDLDETLVHS SFKPISNADF IVPVEIDGTI HQVYVLKRPH VDEFLQRMGQ LFECVLFTAS LAKYADPVAD LLDRWGVFRA RLFRESCVFH RGNYVKDLRS LGRELSKVII VDNPSASYIF HPENAVPVQS WFDDMTDTL LDLIPFFEGL SREDDVYSML HRLCNR
Source	Baculovirus
Target Names	CTDSPL
Protein Names	Recommended name: CTD small phosphatase-like protein Short name= CTDSP-like EC= 3.1.3.16 Alternative name(s): Carboxy-terminal domain RNA polymerase II polypeptide A small phosphatase 3 NIF-like protein Nuclear LIM interactor
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.