



Recombinant Positive regulator of RepFIC repA1 expression (repL)

Product Code	CSB-EP350471SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P62551
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	MPGKVQDFFL CSLLLRIVSA GWCD
Source	E.coli
Target Names	repL
Protein Names	Recommended name: Positive regulator of RepFIC repA1 expression Alternative name(s): RepA1 leader peptide
Expression Region	1-24
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.