



Recombinant Human Nucleoporin-like protein 2 (NUPL2)

Product Code	CSB-MP016211HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15504
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAICQFFLQG RCRFGDRCWN EHPGARGAGG GRQQPQQQPS GNNRRGWNTT SQRYSNVIQP SSFSKSTPWG GSRDQEKPYP SSFD SGASTN RKEGFG LSEN PFASLSPDEQ KDEKKLLEGI VKDMEVWESS GQWMFSVYSP VKKKPNISGF TDISPEELRL EYHNFLTENN LQSYLNSVQR LINQWRNRVN ELKSLNISTK VALLSDVKDG VNQAAPAFGF GSSQAATFMS PGFPVNNSSS DNAQNFSFKT NSGFAAASSG SPAGFGSSPA FGAAASTSSG ISTSAPAFGF GKPEVTS AAS FSFKSPAASS FGSPGFSGLP ASLATGPVRA PVAPAFGGGS SVAGFGSPGS HSHTAFSKPS SDTFGNSSIS TSLSASSSII ATDNVLF TPR DKLTVEELEQ FQSKKFTLGK IPLKPPPLEL LNV
Source	Mammalian cell
Target Names	NUPL2
Protein Names	Recommended name: Nucleoporin-like protein 2 Alternative name(s): NLP-1 NUP42 homolog Nucleoporin hCG1
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.