



# Recombinant Rat Endophilin-A3 (Sh3gl3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP021232RA  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O35180  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSVAGLKKQF HKASQLFSEK ISGAEGTKLD EEFLDMEKKI DITSKAVAEI<br>LSKATEYLQP NPAYRAKLG M LNTMSKLRGQ VKATGYPQTE GLLGDCMLKY<br>GRELGEDSAF GNSLVEVGEA LKLMAEVKDS LDINVKQTFI DPLQLLQDKD<br>LKEIGHHLRK LEGRRLDYDY KKKRVGKIPE EEIRQAVEKF EESKELAERS<br>MFNFLENDVE QVSQLAVFVE AALDYHRQST EILQELQNK L ELRIALASQV<br>PRRDYMPKPV NTSSTNANGV EPSSSSKLTG TDIPSDQPCC RGLYDFEPEN<br>EGELGFKEGD IITLTNQIDE NWYEGMLRGE SGFFPINYVE VIVPLPR |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | Sh3gl3  |
| <b>Protein Names</b>     | Recommended name: Endophilin-A3 Alternative name(s): Endophilin-3 SH3 domain protein 2C SH3 domain-containing GRB2-like protein 3 SH3p13  |
| <b>Expression Region</b> | 1-347   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.       |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |