



Recombinant Pan troglodytes Thyroxine-binding globulin (SERPINA7)

Product Code	CSB-BP351069EQV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61640
Product Type	Recombinant Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Purity	>85% (SDS-PAGE)
Sequence	ASPEGKVTAC HSSQPNATLY KMSSINADFA FNLYRRFTVE TPDKNIFFSP VSISAALVML SFGACCSTQT EIVETLGFNL TDTPMVEIQH GFQHLCISLN FPKKELELQI GNALFIGKHL KPLAKFLNDV KTLYETEVFS TDFSNISAAK QEINSHVEMQ TKGKVVGLIQ DLKPNTIMVL VNYIHFKAQW ANPFDPSKTE DSSSFLIDKT TTVQVPMMHQ MEQYYHLVDM ELNCTVLQMD YSKNALALFV LPKEGQMESV EAAMSSKTLK KWNRLQKGW VDLFVPKFSI SATYDLGATL LKMGIQHAYS ENADFSGLTE DNGLKLSNAA HKAVLHIGEK GTEAAAVPEV ELSDQPENTF LHPIIQIDRS FMLLILERST RSILFLGKVV DPTEA
Source	Baculovirus
Target Names	SERPINA7
Protein Names	Recommended name: Thyroxine-binding globulin Alternative name(s): Serpin A7 T4-binding globulin
Expression Region	21-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.