



Recombinant Lipase

Product Code	CSB-MP351101RDC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61871
Product Type	Recombinant Protein
Immunogen Species	Rhizopus niveus
Purity	≥85% (SDS-PAGE)
Sequence	DDNLV GGMTLDLPSD APPISLSSST NSASDGGKVV AATTAQIQEF TKYAGIAATA YCRSVVPGNK WDCVQCQKWV PDGKIITFT SLLSDTNGYV LRSDKQKTIY LVFRGTNSFR SAITDIVNF SDYKPVKGAK VHAGFLSSYE QVVNDYFPVV QEQLTAHPTY KVIVTGHSLG GAQALLAGMD LYQREPRLSP KNLSIFTVGG PRVGNPTFAY YVESTGIPFQ RTVHKRDIVP HVPPQSFGFL HPGVESWIKS GTSNVQICTS EIETKDCSNS IVPFTSILDH LSYFDINEGS CL
Source	Mammalian cell
Protein Names	Recommended name: Lipase EC= 3.1.1.3 Alternative name(s): Lipase II RNL Triacylglycerol lipase
Expression Region	96-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.