



Recombinant Macrolide-specific efflux protein macA (macA)

Product Code	CSB-EP351391EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P64176
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	PVPTYQTLI VRPGDLQQSV LATGKLDALR KVDVGAQVSG QLKTLSSVAIG DKVKKDQLLG VIDPEQAENQ IKEVEATLME LRAQRQQAEE ELKLARVTYS RQQRLAQTQA VSLQDLDTAA TEMAVKQAQI GTIDAQIKRN QASLDTAKTN LDYTRIVAPM AGEVTQITTL QGQTVIAAQQ APNILTADM STMLVKAQVS EADVIHLKPG QKAWFTVLGD PLTRYEGQIK DVLPTPEKVN DAIFYARFE VPNPNGLLRL DMTAQVHIQL TDVKNVLTIP LSALGDPVGD NRYKVKLLRN GETREREVTI GARNDTDVEI VKGLEAGDEV VIGEAKPGAA Q
Source	E.coli
Target Names	macA
Protein Names	Recommended name: Macrolide-specific efflux protein macA
Expression Region	32-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.