



Recombinant Human Heterogeneous nuclear ribonucleoprotein R (HNRNPR)

Product Code	CSB-EP010614HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43390
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>ANQVNGNAV QLKEEEEPMD TSSVTHTEHY KTLIEAGLPQ KVAERLDEIF QTGLVAYVDL DERAIDALRE FNEEGALSVL QQFKESDLSH VQNKSAFLCG VMKTYRQREK QGSKVQESTK GPDEAKIKAL LERTGYTLDV TTGQRKYGGP PPDVSYSGVQ PGIGTEVFGV KIPRDLYEDE LVPLFEKAGP IWDLRLMMDP LSGQNRGYAF ITFCGKEAAQ EAVKLCDSYE IRPGKHLGVC ISVANNRLFV GSIPKNKTKE NILEEFSKVT EGLVDVILYH QPDDKKKKNRG FCFLEYEDHK SAAQARRRLM SGKVKVWGNV VTVEWADPVE EPDPEVMAKV KVLVFRNLAT TVTEEILEKS FSEFGKLERV KKLKDYAFVH FEDRGAAVKA MDEMNGKEIE GEEIEIVLAK PPKKRKERQ AARQASRSTA YEDYYYHPPP RMPPPIRGRG RGGGRGGYGY PPDYYGYEDY YDDYYGYDYH DYRGGYEDPY YGYDDGYAVR GRGGGRGGRG APPPPRGRGA PPRGRAGYS QRGAPLGPPR GSRGGRGGPA QQQRGRGSRG SRGNRGGNVG GKRKADGYNQ PDSKRRQTNN QQNWGSQPIA QQPLQQGGDY SGNYGYNNDN QEFYQDTYGQ QWK</p>
Source	E.coli
Target Names	HNRNPR
Protein Names	Recommended name: Heterogeneous nuclear ribonucleoprotein R Short name= hnRNP R
Expression Region	2-633
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	<p>This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. This protein has three repeats of quasi-RRM domains that bind to RNAs and also contains a nuclear localization motif. Multiple</p>



alternatively spliced transcript variants have been found for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.