



Recombinant Human Carbohydrate sulfotransferase 10 (CHST10), partial

Product Code	CSB-MP005403HU1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43529
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	TFKDPDVYSAKQEFLFLTTMPEVRKLPPEEKHIPEELKPTGKELPDSQLVQPLVY MERLEL IRNVCRDDALKNLSHTPVSKFVLDRIFVCDKHKILFCQTPKVGNTQWKKVLIVL NGAFSS IEEIPENVVHDHEKNGLPRLSSFSDAEIQKRLKTYFKFFIVRDPFERLISAFKDKF VHNP RFEPWYRHEIAPGIIRKYRRNRRTETRGIQFEDFVRYLGDPNHRWLDLQFGDHII HWVTYV ELCAPCEIMYSVIGHHETLEDDAPYILKEAGIDHLVSYPTIPPGITVYNRTKVEHY FLGI SKRDIRRLYARFEGDFKLFQYQKPDFLLN
Source	Mammalian cell
Target Names	CHST10
Protein Names	Recommended name: Carbohydrate sulfotransferase 10 EC= 2.8.2.-Alternative name(s): HNK-1 sulfotransferase Short name= HNK-1ST Short name= HNK1ST Short name= HuHNK-1ST
Expression Region	28-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.