



# Recombinant Undecaprenyl pyrophosphate synthase (uppS)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP352115SZB   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P60475  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Shigella flexneri   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MMLSATQPLS EKLP AHGCRH VAIIMDGNGR WAKKQGKIRA FGHKAGAKSV<br>RRAVSFAANN GIEALTLYAF SSENWNRPAQ EVSALMELFV WALDSEVKSL<br>HRHNVRLRII GDTSRFNSRL QERIRKSEAL TAGNTGLTLN IAANYGGRWD<br>IVQGVRLAE KVQQGNLQPD QIDEEMLNQH VCMHELAPVD LVIRTGGEHR<br>ISNFWLWQIA YAELYFTDVL WPDFDEQDFE GALNAFANRE RRFGGTEPGD<br>ETA   |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | uppS  |
| <b>Protein Names</b>     | Recommended name: Undecaprenyl pyrophosphate synthase Short name= UPP synthase EC= 2.5.1.31 Alternative name(s): Di-trans,poly-cis-decaprenylcistransferase Ditrans,polycis-undecaprenyl-diphosphate synthase Undecaprenyl diphos   |
| <b>Expression Region</b> | 1-253   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |