



Recombinant *Drosophila hydei* Tubulin beta-2 chain (betaTub85D)

Product Code	CSB-EP352260DLM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61858
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila hydei</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MREIVHIQAG QCGNQIGGKF WEVISDEHCI DATGTYYGDS DLQLERINVV YNEATGAKYV PRAILDLEP GTMDSVRSGA FGQIFRPDNF VFGQSGAGNN WAKGHYTEGA ELVDSVLDVV RKESEGCDCL QGFQLTHSLG GGTGSGMGTL LISKIREEYP DRIMNTFSVV PSPKVSDTVV EPYNATLSVH QLVENTDETY CIDNEALYDI CFRTLKLTTP TYGDLNHLVS ATMSGVTTCL RFPGQLNADL RKLAVNMVPF PRLHFFMPGF APLTSRGSQQ YRALTVPELT QQMFDKNNMM AACDPRHGRY LTVAAIRGR MSMKEVDEQM LNIQKNSSF FVEWIPNNCK TAVCDIPPRG LKMSATFIGN STAIQELFKR VSEQFTAMFR RKAFLHWYTG EGMDMEFTE AESNMNDLVS EYQQYQEATA DEEGEFDEDE EGGGDE
Source	E.coli
Target Names	betaTub85D
Protein Names	Recommended name: Tubulin beta-2 chain Alternative name(s): Beta-2-tubulin
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.