



Recombinant Phospho-2-dehydro-3-deoxyheptonate aldolase, Trp-sensitive (aroH)

Product Code	CSB-EP352609SZB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P59736
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MNRTDELRTA RIESLVTPAE LALRYPVTPG VATHVTDSSRR RIEKILNGED KRLLVIIIGPC SIHDLTAAME YATRLQSLRN QYQSRLEIVM RYFYEKPRTV VGWKGLISDP DLNGSYRVNH GLELARKLLL QVNELGVPTA TEFLDMVTGQ FIADLISWGA IGARTTESQI HREMASALSC PVGFKNGTDG NTRIAVDAIR AARASHMFLS PDKNGQMTIY QTSGNPYGYHI IMRGGKKNPY HADDIAAACD TLHEFDLPEH LVVDFSHGNC QKQHRRQLEV CEDICQQIRN GSTEIAGIMA ESFLREGTQK IVGGQPLTYG QSITDPCLGW EDTERLVEKL AFAVDTRF
Source	E.coli
Target Names	aroH
Protein Names	Recommended name: Phospho-2-dehydro-3-deoxyheptonate aldolase, Trp-sensitive EC= 2.5.1.54 Alternative name(s): 3-deoxy-D-arabino-heptulosonate 7-phosphate synthase DAHP synthase Phospho-2-keto-3-deoxyheptonate aldolase
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.