



Recombinant Mouse Phosphoglycerate mutase 2 (Pgam2)

Product Code	CSB-EP017835MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O70250
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MTTHRLVMVR HGESLWNQEN RFCGWFDAEL SEKGAEAAKR GATAIKDAKI EFDICYTSVL KRAIRTLWTI LDVTDQMWVP VVRTWRLNER HYGGLTGLNK AETAAKHGEE QVKIWRRSFD TPPPMPDEKH NYYTSISKDR RYAGLKPEEL PTCESLKDTI ARALPFWNEE IAPKIKAGQR VLIAAHGNSL RGIVKHLEGM SDQAIMELNL PTGIPIVYEL DQNLKPTKPM RFLGDEETVR KAMEAVAAQG KAK
Source	E.coli
Target Names	Pgam2
Protein Names	Recommended name: Phosphoglycerate mutase 2 EC= 3.1.3.13 EC= 5.4.2.1 EC= 5.4.2.4 Alternative name(s): BPG-dependent PGAM 2 Muscle-specific phosphoglycerate mutase Phosphoglycerate mutase isozyme M Short name= PG
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	Phosphoglycerate mutase (PGAM) catalyzes the reversible reaction of 3- phosphoglycerate (3-PGA) to 2-phosphoglycerate (2-PGA) in the glycolytic pathway. The PGAM is a dimeric enzyme containing, in different tissues, different proportions of a slow-migrating muscle (MM) isozyme, a fast-migrating brain (BB) isozyme, and a hybrid form (MB). This gene encodes muscle-specific PGAM subunit. Mutations in this gene cause muscle phosphoglycerate mutase efficiency, also known as glycogen storage disease X.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.