



Recombinant Rat Plasminogen (Plg)

Product Code	CSB-EP018188RA-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q01177
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	<p>D SLDGYVSTQG ASLHSLTKKQ LAAGSIADCL AKCEGETDFI CRSFQYHSKE QQCVIMAENS KTSSIIRMRD VILFEKRVYL SECKTGIGKG YRGTMSTKTK GVTCQKWSDT SPHVPKYSPS THPSEGLEEN YCRNPDNDEQ GPWCYTTPD QRYEYCNIFE CEEECMYCSG EKYEKISKI MSGLDCQSWD SQSPHAHGYI PAKFPSKNLK MNYCRNPDGE PRPWCFTTDP NKRWEYCDIP RCTTPPPPPG PTYQCLKGRG ENYRGTVSVT ASGKTCQRWS EQTPHRHNRT PENFPCKNLE ENYCRNPDGE TAPWCYTTPD QLRWEYCEIP SCGSSVSPDQ SDSSVLPEQT PVVQECYQGN GKSRYGTSST TNTGKKCQSW VSMTPHSHSK TPANFPDAGL EMNYCRNPDN DQRGPWCFTT DPSVRWEYCN LKRCSETGGG VAESAIVPQV PSAPGTSETD CMYGNGKEYR GKTAVTAAGT PCQEWAQEP HSHRIFTPQT NPRAGLEKNY CRNPDGDVNG PWCYTMNPRK LYDYCNIPLC ASLSSFECGK PQVEPKKCPG RVVGGCVANP HSWPWQISLR TRFSGQHFCG GTLISPEWVL TAAHCLEKSS RPEFYKVLG AHEERILGSD VQQIAVTKLV LEPNDADIAL LKLSRPATIT DNVIPACLPS PNYVVADRTL CYITGWGETK GTPGAGRLKE AQLPVIENTKV CNRAEYLNLR VKSTELCAGH LAGGIDSCQG DSGGPLVCFE KDKYILQGVV SWGLGCARPV KPGVYVRSR YVNWIEREMR ND</p>
Source	E.coli
Target Names	Plg
Protein Names	Recommended name: Plasminogen EC= 3.4.21.7 Cleaved into the following 5 chains: 1. Plasmin heavy chain A 2. Activation peptide 3. Angiostatin 4. Plasmin heavy chain A, short form 5. Plasmin light chain B
Expression Region	20-812
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.