



# Recombinant Human Transcription factor 21 (TCF21)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP023297HU                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                           |
| <b>Uniprot No.</b>       | O43680                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MSTGSLSDVE DLQEVEMLEC DGLKMDSNKE FVTSNESTEE SSNCENGSPQ<br>KGRGGLGKRR KAPTKKSPLS GVSQEGKQVQ RNAANARERA<br>RMRVLSKAFS RLKTTLPWVP PDKLSKLDL LRLASSYIAH LRQILANDKY<br>ENGYIHPVNL TWPFMVAGKP ESDLKEVVTA SRLCGTTAS                                                                                                                                                                                                |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Target Names</b>      | TCF21                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Transcription factor 21 Short name= TCF-21 Alternative name(s): Capsulin Class A basic helix-loop-helix protein 23 Short name= bHLHa23 Epicardin Podocyte-expressed 1 Short name= Pod-1                                                                                                                                                                                                   |
| <b>Expression Region</b> | 1-179                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Details</b>    | TCF21 encodes a transcription factor of the basic helix-loop-helix family. The TCF21 product is mesoderm specific, and expressed in embryonic epicardium, mesenchyme-derived tissues of lung, gut, gonad, and both mesenchymal and glomerular epithelial cells in the kidney. Two transcript variants encoding the same protein have been found for this gene.                                              |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                             |