



# Recombinant Human Striatin (STRN)

<b>Product Code</b>	CSB-MP022880HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O43815
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>MDEQAGPGVF FSNNHGAGG AKGLGPLAEA AAAGDGAAAA  GAARAQYSLP GILHFLQHEW ARFEVERAQW EVERAELQAQ IAFLQGERKG  QENLKKDLVR RIKMLEYALK QERAKYHKLK YGTELNQGDM KPPSYDSDEG  NETEVQPQQN SQLMWKQGRQ LLRQYLQEVG YTDITLDVKS KRVRALLGFS  SDVTDREDDK NQDSVVGTE AEVKETAMIA KSELTDSASV LDNFKFLESA  AADFSEDED DDVDGREKSV IDTSTIVRKK ALPDSGEDRD TKEALKEFDF  LVTSEEGDNE SRSAGDGTDW EKEDQCLMPE AWNVDQGVIT  KLKEQYKKER KGKKGVKRPN RSKLQDMLAN LRDVDELPSL QPSVGSPPSRP  SSSRLPEHEI NRADEVEALT FPPSSGKFSI MGADEALESE LGLGELAGLT  VANEADSLTY DIANNKDALR KTWNPKFTLR SHFDGIRALA FHPIEVLIT  ASEDHTLKMW NLQKTAPAKK STSLDVEPIY TFRAHKGPVL CVVMSSNGEQ  CYSGGTDGLI QGWNTTNPNI DPYDSYDPSV LRGPLLGHTD AVWGLAYSAA  HQRLLSCSAD GTLRLWNTTE VAPALSVFND TKELGIPASV DLVSSDP SHM  VASFSGYTS IFNMQQR I LTLESNVDTT ANSSCQINRV ISHPTLPISI  TAHEDRHIKF YDNNTGKLIH SMVAHLEAVT SLAVDPNGLY LMSGSHDCSI  RLWNLESKTC IQEFTAHRKK FEESIHDVAF HPSKCYIASA GADALAKVVF</p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	STRN
<b>Protein Names</b>	Recommended name: Striatin
<b>Expression Region</b>	1-780
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.