



Recombinant Response regulator uvrY (uvrY)

Product Code	CSB-EP353627EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P66797
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MINVLLVDDH ELVRAGIRRI LEDIKGIKVV GEASCGEDAV KWCRANAVDV VLMDMSMPGI GGLEATRKIA RSTADVKIIM LTVHTENPLP AKVMQAGAAG YLSKGAAPQE VVSAIRSVYS GQRYIASDIA QQMALSQIEP EKTESPFASL SERELQIMLM ITKGQKVNEI SEQLNLSPKT VNSYRYRMFS KLNHGDVEL THLAIRHGLC NAETLSSQ
Source	E.coli
Target Names	uvrY
Protein Names	Recommended name: Response regulator uvrY
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.