



Recombinant Escherichia coli Uncharacterized protein YhcM (yhcM)

Product Code	CSB-EP353749ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P64612
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MQSVTPTSQY LKALNEGSHQ PDDVQKEAVS RLEIIYQELI NSTPPAPRTS GLMARVGKLW GKREDTKHTP VRGLYMWGGV GRGKTWLM DL FYQSLPGERK QRLHFHRFML RVHEELTALQ GQTDPLEIIA DRFKAETDVL CFDEFFVSDI TDAMLLGGLM KALFARGITL VATSNIPPDE LYRNLQRAR FLPAIDAIKQ HCDVMNVDAG VDYRLRRTLQ AHLWLSPLHD ETRAQMDKWLW LALAGGKREN SPTLEINHRP LATMGVENQT LAVSFTTLCV DARSQHDYIA LSRLFHTVML FDVPVMTRLM ESEARRFIAL VDEFYERHVK LVSVAEVP LY EIQGDRLKF EFQRCL SRLQ EMQSEEYLKR EHLAG
Source	E.coli
Target Names	zapE
Protein Names	Recommended name: Uncharacterized protein YhcM
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.