



Recombinant Human Leucine carboxyl methyltransferase 2 (LCMT2)

Product Code	CSB-MP012809HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60294
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MGPRSRERRA GAVQNTNDSS ALSKRSLAAR GYVQDPFAAL LVPGAARRAP LIHRGYVRA RAVRHCVRAF LEQIGAPQAA LRAQILSLGA GFDSLYFRLK TAGRLARA AV WEVDFPDVAR RKAERIGETP ELCALTGPFE RGE PASALCF ESADYCILGL DLRQLQRVEE ALGAAGLDAA SPTLLLAEAV LTYLEPESAA ALIAWAAQRF PNALFVVEEQ MRPQDAFGQF MLQHFRQLNS PLHGLERFPD VEAQRRRFLQ AGWTACGA VD MNEFYHCFLP AEERRRVENI EPFDEFEEWH LKCAHYFILA ASRGDTLSHT LVFPSSEAFP RVNPASPSGV FPASVVSSEG QVPNLKRYGH ASVFLSPDVI LSAGGFGEQE GRHCRVSQFH LLSRDCDSEW KGSQIGSCGT GVQWDGRLYH TMLRLSESRV LVLGGRLSPV SPALGVLQLH FFKSEDNTE DLKVTITKAG RKDDSTLCCW RHSTTEVSCQ NQEYLFVYGG RSVVEPVLSD WHFLHVG TMA WVRIPVEGEV PEARHSHSAC TWQGGALIAG GLGASEEPLN SVLFLRPISC GFLWESVDIQ PPITPRYSHT AHVLNGKLLL VGGIWIHSSS FPGVTVINLT TGLSSEYQID TTYVPWPLML HNHTSILLPE EQQLLLLGGG GNCFSFGTYF NPHTVTL DLS SLSAGQ</p>
Source	Mammalian cell
Target Names	LCMT2
Protein Names	<p>Recommended name: Leucine carboxyl methyltransferase 2 EC= 2.1.1.- Alternative name(s): p21WAF1/CIP1 promoter-interacting protein tRNA wybutosine-synthesizing protein 4 homolog Short name= tRNA yW-synthesizing protein 4 homolog</p>
Expression Region	1-686
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The protein encoded by this intronless gene belongs to the highly variable methyltransferase superfamily. This gene is the inferred homolog of the <i>Saccharomyces cerevisiae</i> carboxymethyltransferase gene PPM2 that is essential for the synthesis of the hypermodified guanosine Wybutosine (yW).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.