



Recombinant Human PRAME family member 10 (PRAMEF10)

Product Code	CSB-MP018605HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60809
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLQAPSRLI ELAQSLLRN QFLTIFTLDE LPREVFPLMF MEAFSMRRFE ALKLMVQAWP FLSLPLGSLM KTLHLETQLA VLRGLDTLVA QKVRPRRWKL QVLDLRDVDE NFWTIWVGAR VLSCSPEAMS KRQTVEDCPR MGEHQPLKVF IDLCLKESTL DECLSYLFGW IHYRRGLVHL CCSKVQNYSM PTSSFRNLLE RIYPDSIQEL EVWKKCSLNK TGKFAPYLSQ MSNRELFLA FGYERELYVS VQWPCIPDL SPFLCLYYPQ MLYIKKISNI KEHLEHLLRY LKNPLGAFIF SDAYLTDRDM ECLSQYPSLS QLKELRLIHI LMWTTTHLEPL GVLLKVAAT LKTLVLKDCR IQDPQLRVLL PALSHCSQLT TFNFHGNETS MNALKDLLRH TGGLSKLGL LYPAPLES LD YKGHVNWEIL TPIRAELMRT LREVRQPKRI FFGPVPCPNC GSWPSEKVDV HLCS
Source	Mammalian cell
Target Names	PRAMEF10
Protein Names	Recommended name: PRAME family member 10
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.