



# Recombinant Human Coiled-coil domain-containing protein 22 (CCDC22)

<b>Product Code</b>	CSB-EP004681HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O60826
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEEADRILIH SLRQAGTAVP PDVQTLRAFT TELVVEAVVR CLRVINPAVG SGLSPLLPLA MSARFRLAMS LAQACMDLGY PLELGYQNFL YPSEPDRLDL LLFLAERLPT DASEDADQPA GDSAILLRAI GSQIRDQLAL PWVPPHLRTP KLQHLQGSAL QKPFHASRLV VPELSSRGEP REFQASPLLL PVPTQVPQPV GRVASLLEHH ALQLCQQTGR DRPGDEDWVH RTSRLPPQED TRAQRQRLQK QLTEHLRQSW GLLGAPIQAR DLGELLQAWG AGAKTGAPKG SRFTHSEKFT FHLEPQAQAT QVSDVPATSR RPEQVTWAAQ EQELESREQ LEGVNRSIEE VEADMKTLGV SFVQAESECR HSKLSTAERE QALRLKSRV ELLPDGTANL AKLQLVVENS AQRVIHLAGQ WEKHRVPLLA EYRHLRKLQD CRELESSRRL AEIQELHQSVA RAAAEARRK EEVYKQLMSE LETLPRDVSR LAYTQRILEI VGNIRKQKEE ITKILSDTKE LQKEINSLSG KLDRTFAVTD ELVFKDAKKD DAVRKAYKYL AALHENCSQL IQTIEDTGTI MREVRDLEEQ IETELGKCTL SNLEKIREDY RALRQENAGL LGRVREA
<b>Source</b>	E.coli
<b>Target Names</b>	CCDC22
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 22
<b>Expression Region</b>	1-627
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.