



Recombinant Human Keratocan (KERA)

Product Code	CSB-BP012149HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60938
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	RSVRQVYEVH DSDDWTIHDF ECPMECF CPP SFPTALYCEN RGLKEIPAIP SRIWYLYLQN NLIETIPEKP FENATQLRWI NLNKNKITNY GIEKGALSQ KLLFLFLED NELEEVPSPL PRSLEQLQLA RNKVSRI PQG TFSNLENLTL LDLQNNKLVD NAFQRDTFKG LKNLMQLNMA KNALRNMP PR LPANTMQLFL DNNSIEGIPE NYFNVIPKVA FLRLNHNKLS DEGLPSRGFD VSSILDQLS HNQLTKVPRI SAHLQHLHLD HNKIKSVNVS VICPSPMLP AERDSFSYGP HLRYLRLDGN EIKPPIPMAL MTCFRLLQAV II
Source	Baculovirus
Target Names	KERA
Protein Names	Recommended name: Keratocan Short name= KTN Alternative name(s): Keratan sulfate proteoglycan keratocan
Expression Region	21-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.