



Recombinant Mouse Midline-1 (Mid1)

Product Code	CSB-MP013819MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O70583
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	METLESELTCPICLELFEDP LLLPCAHS LC FNCAHRILVS HCATNEPVES INAFQCPTCR HVITLSQRGL DGLKRNVT LQ NIIDRFQKAS VSGPNSPSET RRERAFDANT MSSAEKVLCQ FCDQDPAQDA VKTCVTCEVS YCDECLKATH PNKKPFTGHR LIEPIPSHI RGLMCLEHED EKVNMVCVTD DQLICALCKL VGRHRDHQVA ALSERYDKLK QNLESNLTNL IKRNTELETL LAKLIQTCQH VEVNASRQEA KLTEECDLLI EIIQRRRQII GTKIKEGKVI RLRKLAQQIA NCKQCLERSA SLISQAEHSL KENDHARFLQ TAKNITERVS MATASSQVLI PEINLNDTFD TFALDFSREK KILLECLDYLT APNPPAIREE LCTASYDTIT VHWTSEDEFS VVS YELQYTI FTGQANVVNV ACDGTCLLGS AGLCNSADSW MIVPNIKQNH YTVHGLQSGT KYIFTVKAIN QAGSRSSPEPG KLKTNSQPFR LDPKSAHRKL KVSHDNLTVE RDESSSKKSH APERFAGQGS YGVAGNVFID SGRHYWEVVT SGSTWYAIGL AYRSAPKHEW IGKNAASWAL CRCHNHWA VR HDGKETPIAP APHLRRVGV L LDYDNGSIAF YDALSSVHLH TFHAALAQPV CPTFTVWNKC LTIVTGLPIP DHLDCTEQRP
Source	Mammalian cell
Target Names	Mid1
Protein Names	Recommended name: Midline-1 EC= 6.3.2.- Alternative name(s): Midin Tripartite motif-containing protein 18
Expression Region	1-680
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is a member of the tripartite motif (TRIM) family, also known as the RING-B box-coiled coil (RBCC) subgroup of RING finger proteins. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein forms homodimers which associate with microtubules in the cytoplasm. The protein is likely involved in the formation of multiprotein structures acting as anchor points to microtubules. Mutations in this gene have been associated with the X-linked form of Opitz syndrome, which is characterized by midline abnormalities such as cleft lip, laryngeal cleft, heart defects, hypospadias, and agenesis of the corpus callosum. This gene was also the first example of a gene subject to X inactivation in human while escaping it in mouse. Several different transcript variants are generated by alternate



splicing; however, the full-length nature of some of the variants has not been determined.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.