



# Recombinant UDP-N-acetylmuramate--L-alanine ligase (murC)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP354543EGX  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | P65470   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MNTQQLAKLR SIVPEMRRVR HIHFVGIGGA GMGGIAEVLA NEGYQISGSD<br>LAPNPVTQQL MNLGATIYFN HRPENVRDAS VVVVSSAISA DNPEIVAHE<br>ARIPVIRRAE MLAELMRFRH GIAIAGTHGK TTTTAMVSSI YAEAGLDPTF<br>VNGGLVKAAG VHARLGHGRY LIAEADESDA SFLHLQPMVA IVTNIADHM<br>DTYQGDFENL KQTFINFLHN LPFYGRAVMC VDDPVIRELL PRVGRQTTTY<br>GFSEDADVRV EDYQQIGPQG HFTLLRQDKE PMRVTLNAPG RHNALNAAAA<br>VAVATEEGID DEAILRALES FQGTGRRFDF LGEFPLEPVN GKSGTAMLVD<br>DYGHHPTVD ATIKAAARAGW PDKNLVMLFQ PHRFTRTRDL YDDFANVLTQ<br>VDTLLMLEVY PAGEAPIGA DSRSLCRTIR GRGKIDPILV PDPAQVAEML<br>APVLTGNDLI LVQGAGNIGK IARSLAEIKL KPQTPEEEQH D |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | murC   |
| <b>Protein Names</b>     | Recommended name: UDP-N-acetylmuramate--L-alanine ligase EC= 6.3.2.8<br>Alternative name(s): UDP-N-acetylmuramoyl-L-alanine synthetase   |
| <b>Expression Region</b> | 1-491  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |