



Recombinant Ribosomal large subunit pseudouridine synthase D (rluD)

Product Code	CSB-MP354612SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P65835
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	AQRVQLTAT VSENQLGQRL DQALAEMFPD YSRSRKEWI LDQRVLVNGK VCDKPKEKVL GGEQVAINAE IEEEARFEPQ DIPLDIVYED EDIIVINKPR DLVVHPGAGN PDGTVLNALL HYYPPIADVP RAGIVHRLDK DTTGLMVVAK TVPAQTRLVE SLQRREITRE YEAVAIGHMT AGGTVDEPIS RHPTKRTHMA VHPMGKPAVT HYRIMEHFRV HTRLRLRLET GRTHQIRVHM AHITHPLVGD PVYGGRRPPP KGASEAFIST LRFKDRQALH ATMLRLYHPI SGIEMEWHP IPQDMVELIE VMRADFEEHK DEVDWL
Source	Mammalian cell
Target Names	rluD
Protein Names	Recommended name: Ribosomal large subunit pseudouridine synthase D EC=5.4.99.23 Alternative name(s): 23S rRNA pseudouridine(1911/1915/1917) synthase rRNA pseudouridylylate synthase D rRNA-uridine isomerase D
Expression Region	2-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.