



# Recombinant Human Cold shock domain-containing protein E1 (CSDE1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP006040HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O75534  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSFDPNLLHN NGHNGYPNGT SAALRETGVI EKLLTSYGFI QCSEQRARLF<br>FHCSQYNGNL QDLKVGDDVE FEVSSDRRTG KPIAVKLVKI KQEILPEERM<br>NGQVVCAPPH NLESKSPAAP GQSPTGVCY ERNGEVFYLT YTPEDVEGNV<br>QLETGDKINF VIDNNKHTGA VSARNIMLLK KKQARCQGVV CAMKEAFGFI<br>ERGDVVKEIF FHYSEFKGDL ETLQPGDDVE FTIKDRNGKE VATDVRLLPQ<br>GTVIFEDISI EHFEGTVTKV IPKVPSKNQN DPLPGRIKVD FVIPKELPFG<br>DKDTKSKVTL LEGDHVRFNI STDRRDKLER ATNIEVLSNT FQFTNEAREM<br>GVIAAMRDGF GFIKCVDRDV RMFFHFSEIL DGNQLHIAD EFTVVPDML<br>SAQRNHAIRI KKLPGKTVSF HSHSDHRFLG TVEKEATFSN PKTTSPNKGK<br>EKEAEDGIIA YDDCGVKLTI AFQAKDVEGS TSPQIGDKVE FSISDKQRPG<br>QQVATCVRLG GRNSNSKRLG GYVATLKDNF GFJETANHDK EIFFHYSEFS<br>GDVDSLELGD MVEYSLSKGG GNVSAEKVN KTHSVNGITE EADPTIYSGK<br>VIRPLRSVDP TQTEYQGMIE IVEEGDMKGE VYPPFGVGMMA NKGDCCLQKGE<br>SVKFQLCVLG QNAQTMAYNI TPLRRATVEC VKDQFGFINY EVGDSKLLFF<br>HVKEVQDGIE LQAGDEVEFS VILNQRGK SACNVWRVCE GPKAVAAPRP<br>DRLVNRLKNI TLDDASAPRL MVLRRQPRGPD NSMGFGAERK IRQAGVID |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | CSDE1   |
| <b>Protein Names</b>     | Recommended name: Cold shock domain-containing protein E1 Alternative name(s): N-ras upstream gene protein Protein UNR  |
| <b>Expression Region</b> | 1-798   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,   |



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.