



Recombinant Human TIP41-like protein (TIPRL)

Product Code	CSB-EP023569HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75663
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMIHGFQSSH RDFCFGPWKL TASKTHIMKS ADVEKLADEL HMPSLPEMMF GDNVLRIQHG SGFGIEFNAT DALRCVNNYQ GMLKVACAAE WQESRTEGEH SKEVIKPYDW TYTTDYKGT L GESLKLKVV PTTDHIDTEK LKAREQIKFF EEVLLFEDEL HDHGVSSLSV KIRVMPSSFF LLLRFFLRID GVLIRMNDTR LYHEADKTYM LREYTSRESK ISSLMHVPPS LFTEPNEISQ YLPIKEAVCE KLIFPERIDP NPADSQKSTQ VE
Source	E.coli
Target Names	TIPRL
Protein Names	Recommended name: TIP41-like protein Alternative name(s): Putative MAPK-activating protein PM10 Type 2A-interacting protein Short name= TIP
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.