




Recombinant Human Renin receptor (ATP6AP2)

Product Code	CSB-YP002384HU
Relevance	Functions as a renin and prorenin cellular receptor. May mediate renin-dependent cellular responses by activating ERK1 and ERK2. By increasing the catalytic efficiency of renin in AGT/angiotensinogen conversion to angiotensin I, it may also play a role in the renin-angiotensin system (RAS).
Abbreviation	Recombinant Human ATP6AP2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75787
Alias	ATPase H(+)-transporting lysosomal accessory protein 2 ATPase H(+)-transporting lysosomal-interacting protein 2 ER-localized type I transmembrane adaptor Embryonic liver differentiation factor 10 N14F Renin/prorenin receptor Vacuolar ATP synthase membrane sector-associated protein M8-9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	NEFSILKSPGSSVFRNGNWPIPIGERIPDVAALSMGFVSKEDLSWPGLAVGNLF HRPRATVMVMVKGVNKLALPPGSVISYPLENAVPFSLDSVANSIHSLFSEETPV VLQLAPSEERVYMGKANSVFEDLSVTLRQLRNRLFQENSVLSSLPLNLSLRN NEVDLLFLSELQVLHDISSLLSRHKHLAKDHSPDLYSLELAGLDEIGKRYGEDS EQFRDASKILVDALQKFADDMYSLYGGNAVVELVTVKSFDTSLIRKTRTILEAK QAKNPASPYNLAYKYNFEYSVVFNMVLWIMIALALAVIITSYNIWNMDPGYDSII YRMTNQKIRMD
Research Area	Signal Transduction
Source	Yeast
Target Names	ATP6AP2
Protein Names	Recommended name: Renin receptor Alternative name(s): ATPase H(+)-transporting lysosomal accessory protein 2 ATPase H(+)-transporting lysosomal-interacting protein 2 ER-localized type I transmembrane adaptor Embryonic liver differentiation f
Expression Region	17-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	40.0kDa

**Protein Length**

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.